

# PM-600

DIGITAL PD METER



## PM-600 Specifications

Measurement range	47 mm to 83 mm
Measurement step	0.5 mm
Display method	LCD
Viewing point distance	30 cm, 35 cm, 40 cm, 50 cm, 65 cm, 1 m, 2 m, infinity
Auto power off	Turned off after approximately 3 minutes
Single-eye shutter	LCD method
Hairline	LCD method
Power supply	4.5 V (SUM-3 x3)
Dimensions & Weight	146 (W) x 225 (H) x 50 (D) mm / 590 g 5.7 (W) x 8.9 (H) x 2.0 (D)" / 1.3 lbs.

## Key Features

- Instant PD calculation of each viewing point distance
- Power saving function: automatic power off



Eye & Health Care  
**NIDEK CO., LTD.**

### HEAD OFFICE

34-14 Maehama, Hiroishi  
Gamagori, Aichi 443-0038, Japan  
Telephone : 81-533-67-6611  
Facsimile : 81-533-67-6610  
URL : <http://www.nidek.co.jp>

### TOKYO OFFICE

(International Div.)  
3F Sumitomo Fudosan Hongo Bldg.,  
3-22-5 Hongo, Bunkyo-ku, Tokyo,  
113-0033 Japan  
Telephone : 81-3-5844-2641  
Facsimile : 81-3-5844-2642  
URL : <http://www.nidek.com>

### NIDEK INC.

47651 Westinghouse Drive  
Fremont, CA 94539, U.S.A.  
Telephone : 1-510-226-5700  
: 1-800-223-9044 (US only)  
Facsimile : 1-510-226-5750  
URL : <http://www.usa.nidek.com>

### NIDEK SOCIÉTÉ ANONYME

Europarc  
13, rue Auguste Perret  
94042 Créteil, France  
Telephone : 33-1-49 80 97 97  
Facsimile : 33-1-49 80 32 08  
URL : <http://www.nidek.fr>

### NIDEK TECHNOLOGIES SRL

Via dell'Artigianato, 6 / A  
35020 Albignasego (Padova), Italy  
Telephone : 39 049 8629200 / 8626399  
Facsimile : 39 049 8626824  
URL : <http://www.nidektechnologies.it>



\*Specifications and design are subject to change without notice for improvement.



Printed on environment-friendly recycled paper.