



# CB for Windows

for RT-6100/TS-610

Stay connected with the patients from a remote location.\*1  
NIDEK CB for Windows will improve the flexibility and efficiency of your operation.



*Expanding possibilities,*

*Enhancing refractions*




CB for Windows is an optional control software for the Intelligent Refractor RT-6100 and the Tabletop Refraction System TS-610. This software delivers the comprehensive examination capabilities of the existing control console through a Windows computer.

Experience the advantage of new possibilities provided by NIDEK.



## Control software for RT-6100/TS-610

Featuring a highly visible measurement screen and digital keypad on the computer display, the lens exchanges are as seamless as using a control console. Practical functions are available including pre-programmed refraction settings for consistent and reliable examinations.



## More technology, More efficiency

CB for Windows enables remote examinations by using functions and software<sup>\*1</sup> through a Windows computer to control the equipment remotely. Combined with the Fully Assisted Refraction System (FARS), an optional kit for the NIDEK Tabletop Refraction System TS-610, a patient checks with subjective refraction to achieve the full correction values in a facility then the eye care professional provides the final check/adjustment remotely.

<sup>\*1</sup> Operation from a different location can be performed by remote operation of the computer. NIDEK does not include or provide the functions or software necessary for remote communications. The working environment and the specifications/performance of individual computers may impact software usability.

### Selectable user interface

The operator can choose between using the standard control console or a computer/tablet by using CB for Windows.<sup>\*2</sup> Explore the possibilities.



<sup>\*2</sup> Control console and Windows computer/tablet cannot be used/connected simultaneously.

### CB for Windows Specifications

Software	RT-6100 CB for Windows (installation CD)
Standard accessories	License key
Compatible devices	RT-6100, TS-610
Optional accessories for TS-610	Control box tray, Tablet attachment

### CB for Windows System Requirements

Operating system	Windows 10 (32-bit/64-bit) or Windows 11 (64-bit)
Display	1,920 x 1,080 pixels or higher
Disk device	CD-ROM drive (only for installation)
Other	USB2.0: 1 port or more LAN (100BASE-T or higher): 1 port or more Mouse SD card reader (if necessary)



Intelligent Refractor RT-6100



Tabletop Refraction System TS-610



Control box tray for TS-610



Tablet attachment for TS-610

Product/model name: REFRACTOR RT-6100

CB for Windows is optional for the Intelligent Refractor RT-6100 and the Tabletop Refraction System TS-610. CB for Windows does not work without these devices.

Brochure and listed features of the product are intended for non-US practitioners.

Specifications may vary depending on circumstances in each country.

Specifications and design are subject to change without notice.

The computer and tablet described in this brochure are not provided by NIDEK.

Microsoft and Windows are registered trademarks or trademarks of Microsoft Corporation in the United States and other countries.

All LCD images are simulated.

